Quirem Medical B.V.
Zutphenseweg 55
7418 AH Deventer
The Netherlands
www.QuiremMedical.com

1. DESCRIPTION

QuiremSpheres® consists of biocompatible poly-L-lactic acid (PLLA) microspheres containing Holmium-166. The microspheres have a mean diameter of 30 micrometer (97% between 15 to 60 micrometer). Holmium-166 is a high-energy beta emitting isotope for therapeutic use. The maximum energy of the beta particles is 1.85 MeV (50.0%) and 1.77 MeV (48.7%). The maximum range of the emitted beta particles in tissue is 8.7 mm with a mean of 2.5 mm. In addition, Holmium-166 emits primary gamma photons (81 KeV). The half-life is 26.8 hours, which means more than 90% of the radiation is delivered within the first 4 days following the implantation procedure. At planned moment of treatment, the activity per microsphere is approximately 240-375 Bq. QuiremSpheres® is a permanent implant.

QuiremSpheres® is supplied as a patient specific dose in a single capped V-vial. Details on the activity at the reference time and expiry information are indicated on the package labels. Each V-vial comes in a lead container, packed in a type A transport box. The vial and its content should be stored inside its transportation container at room temperature (15-25°C, 59-77°F) until use.

Caution: Upon receiving QuiremSpheres®, check for damage of the packaging. Do not use the product if the V-vial or package is damaged. QuiremSpheres® is for single patient use only. Do not re-sterilize the product and do not use the product after the expiry date and time, as indicated on the label.

Before using QuiremSpheres®, verify the activity of the received vials by using a calibrated well-chamber or dosimeter.

QuiremSpheres® is administered into the hepatic artery via a catheter. QuiremSpheres® will distribute non-uniformly in the liver. This is primarily due to differences in the hepatic arterial flow to the tumour(s) and non-tumour liver tissue, the ratio of vascularity between tumour and non-tumour liver tissue, and tumour size.

QuiremSpheres® microspheres can be visualized in-vivo with MRI and/or SPECT.

2. INTENDED USE

QuiremSpheres® is intended for implantation into hepatic tumours by delivery via the hepatic artery.

3. INDICATIONS FOR USE

QuiremSpheres® is indicated for the treatment of unresectable liver tumours.

4. ACCESSORIES

Accessories that may be used for the implant procedure include

- QuiremSpheres® Delivery Set (QS-D001)
- QuiremSpheres® Customer Kit (QS-C001)

5. INSTRUCTIONS FOR USE

5.1 PATIENT SELECTION AND PRE-TREATMENT TESTING

Patients with non-resectable liver tumours may be considered for treatment with QuiremSpheres®. Patient selection for treatment with QuiremSpheres® requires a medical opinion that control of tumour within the liver will result in patient benefit and that the patient has received all available surgical treatments.

Patient Tests Before Treatment with QuiremSpheres®

The following tests are recommended to be performed before treatment:

- A hepatic angiogram to establish arterial anatomy of the liver.
- A nuclear medicine break-through scan (Intrahepatic Technetium-99m MAA Scan or other surrogate marker) to determine the lung shunt. This is performed through the hepatic artery catheter.
- Serologic tests of liver function to determine the extent of liver function damage.

5.2 RADIATION DOSIMETRY

The recommended average whole-liver absorbed dose is 60 Gy. Due to non-uniform distribution of QuiremSpheres® to the tumour and non-tumour liver tissue, a proportionally higher radiation dose will be delivered to the tumour tissue.

5.3 CALCULATION OF INDIVIDUAL DOSE

The amount of Holmium-166 radioactivity (A) in [MBq] to be administered to a patient to deliver a liver dose LD in [Gy] is calculated using the following formula:

$$A[\text{MBq}] = LD[\text{Gy}] \times LW[\text{kg}] \times 63 [\text{MBq/J}]$$

where LW is the liver weight in [kg] that will be treated. The liver weight is determined based on CT or MR images.

Using the recommended average whole-liver absorbed dose of 60 Gy ($LD = 60$ Gy), the above formula can be simplified into:

$$A[\text{MBq}] = 3781 [\text{MBq/kg}] \times LW[\text{kg}].$$

If this total activity (A) is planned to be delivered to distinct segments in separate administrations, the required activity per segment (A_i) should be calculated according to the following formula:

$$A_i[\text{MBq}] = W_i [\text{kg}]/LW[\text{kg}] \times A[\text{MBq}]$$

where W_i is the weight of the specific segment determined based on CT or MR images. LW and (A) are the previously mentioned liver weight and total required activity.

The sum of the activity of all segments should always be equal to the total activity to be administered to the patient, (A).

5.4 TECHNIQUE FOR MEASURING LUNG SHUNT

When using a Technetium-99m MAA (Macro-Aggregated Albumin) scan to assess lung shunt, the following considerations should be made. Technetium-99m MAA with an activity of 150 MBq (4 mCi) is injected via a catheter into the hepatic artery. The distribution of Technetium-99m MAA is assessed by SPECT/CT or planar scintigraphy. The formula below is used to determine the

fraction that passes through the liver into the lungs (lung shunt fraction). In addition, any extrahepatic deposition should be evaluated, as well as the intrahepatic distribution of Technetium-99m-MAA between tumour and non-tumour liver tissue.

To optimize the predictive value of the Technetium-99m-MAA scan, the catheter tip during Technetium-99m-MAA injection is recommended to be placed at the exact same anatomical position from which QuiremSpheres® will be administered.

Regions of interest are drawn around the entire lung area and the entire liver area. The percentage of lung shunting (L) is calculated using the following formula:

$$L = \frac{\text{counts of total lung}}{\text{counts of total lung} + \text{counts of liver}} \times 100$$

Prediction of Holmium-166 lung absorbed dose

The Holmium-166 lung absorbed dose can be predicted/estimated based on the Technetium-99m MAA lung shunt (L), the amount of Holmium-166 activity (A) to be administered and the mass of the lungs (M_{lung})

$$D_{lung} \text{ [Gy]} = \frac{A[\text{MBq}] \times L/100}{M_{lung}[\text{kg}] \times 63 [\text{MBq/l}]}$$

Lung tissue mass (M_{lung}) can be calculated using a patient-specific lung tissue density and lung volume based on a CT scan, or can be estimated as 1 kg.

Lung shunt should not lead to a lung dose exceeding 30 Gy in a single treatment.

Table 1 – Patient eligibility according to the estimated average lung absorbed dose

Predicted lung dose (based on Tech-99m-MAA scan)	Recommendation
>30 Gy	Do not use QuiremSpheres®

5.5 PREPARING AND EXECUTING THE IMPLANTATION PROCEDURE

It is strongly recommended to use the following accessories for the QuiremSpheres® administration procedure:

- QuiremSpheres® Delivery Set (QS-D001)
- QuiremSpheres® Customer Kit (QS-C001)

For an extended and illustrated description of the preparation and execution of the QuiremSpheres® administration procedure, the reader is referred to the instructions for use of the QuiremSpheres® Customer Kit and Delivery Set (document number MAN-1101-16-01). These instructions should be read and understood in their entirety prior to use.

Transarterial Implantation

For transarterial implantation of QuiremSpheres® a catheter is inserted either via the femoral or the radial artery under x-ray guidance. This should only be performed by a trained interventional radiologist.

Routinely check the catheter tip to ensure it remains in the planned position throughout the implant procedure. The transarterial catheter can be inserted into hepatic artery branches to reduce the risk of reflux of QuiremSpheres® into small arteries supplying, for example, the gut.

Once the catheter has been correctly positioned and is connected to the QuiremSpheres® Delivery Set, the procedure for delivering QuiremSpheres® is as follows:

- The hepatic artery catheter should be placed into the arterial supply of the liver so that the target area of the liver is adequately perfused.
- Small arteries that pass from the common hepatic artery (and sometimes even from the right or left hepatic arteries) to the stomach and duodenum can be coiled to avoid extrahepatic deposition into the stomach and duodenum.
- It is advised to use a catheter with an inner diameter of at least 0.8 mm. If a catheter with smaller diameter is used, occlusion may occur during the delivery of QuiremSpheres®.
- Delivery of QuiremSpheres® must be performed **slowly (< 5 ml per minute)** using saline (0.9%) into the hepatic artery to prevent the microspheres from refluxing back into the hepatic artery and lodging inside the pancreas, stomach and/or other organs. During the delivery procedure, the catheter should be flushed with saline (0.9%) at regular intervals to prevent blockage.

Caution: QuiremSpheres® must be delivered slowly at a rate of no more than 5 ml per minute. Rapid delivery may cause reflux back down the hepatic artery and into other organs.

Caution: During administration, stasis of blood flow and potential microspheres reflux should be checked intermittently by administration of contrast agent to prevent inadvertent delivery of microspheres.

Once the administration procedure is completed, the catheter can be removed.

Radiological Placement of Catheter

The radiologist must be familiar with the frequent arterial abnormalities in the blood supply to the liver and from the liver to the gut. Every attempt should be made to deliver the QuiremSpheres® microspheres into the hepatic arteries in such a way that radiation is administered to the target area in the liver only. If the tumours are limited to one lobe, the catheter can be inserted selectively into the lobar artery supplying the target lobe, thereby sparing the healthy lobe.

It is essential that QuiremSpheres® is not delivered to other organs, in particular the pancreas, stomach or duodenum. In order to prevent the delivery of QuiremSpheres® to duodenum and stomach, the catheter must be placed well distal to the gastro-duodenal artery (GDA) and any other artery that is supplying blood to the gut. If there is any possibility of QuiremSpheres® passing down the GDA, the implantation must not proceed. It may be preferable to block the GDA with an intraluminal coil or other method to prevent QuiremSpheres® from flowing to the duodenum.

Note: Virtually all complications from QuiremSpheres® arise from the inadvertent delivery of QuiremSpheres® into small blood vessels that go to the pancreas, stomach or duodenum, or inadvertent delivery in the liver itself (i.e. high non-tumour liver absorbed dose).

It is recommended to perform a post-treatment SPECT scan of the upper abdomen to review the intrahepatic distribution and extrahepatic deposition of the administered microspheres.

6. CONTRAINDICATIONS

QuiremSpheres® microspheres are contraindicated in patients who have:

- had previous external beam radiation therapy to the liver;
- ascites or are in clinical liver failure;
- inadequate liver function
- significantly abnormal synthetic and excretory liver function tests (LFTs);
- a lung dose exceeding 30 Gy in a single treatment.
- a pre-assessment angiogram that demonstrates abnormal vascular anatomy that would result in significant reflux of hepatic arterial blood to the stomach, pancreas or bowel;
- been treated with capecitabine within two months prior to treatment, or who will be treated with capecitabine at any time following treatment with QuiremSpheres®;
- complete main portal vein thrombosis;
- inadequate kidney function
- uncorrectable extrahepatic deposition. Activity in the falciform ligament, portal lymph nodes and gallbladder is accepted.

7. WARNINGS

1. Inadvertent delivery of QuiremSpheres® to the gastrointestinal tract or pancreas will cause acute abdominal pain, acute pancreatitis or peptic ulceration.
2. High levels of implanted radiation and/or excessive shunting to the lungs (> 30 Gy in a single session) may lead to radiation pneumonitis.
3. Excessive radiation to the normal liver parenchyma may result in Radioembolization-induced liver disease (REILD).
4. Inadvertent delivery of QuiremSpheres® to the gall bladder may result in cholecystitis.
5. The patient must recover from any major surgical operations before being treated with QuiremSpheres®.

8. PRECAUTIONS

- Only clinical staff that are trained under the QuiremSpheres® training program may order, handle and/or implant QuiremSpheres® microspheres.
- Safety and effectiveness of this device in pregnant women, nursing mothers or children have not been established.
- This product is radioactive. Local regulations must be followed when handling this device.
- Some patients may develop gastritis following treatment. Gastric acid blocking drugs may be used the day before implantation of QuiremSpheres® and continued as needed to reduce gastric complications.
- Patients may experience abdominal pain immediately after administration of QuiremSpheres® and pain relief may be required.

9. ADVERSE EVENTS

The common adverse events after receiving QuiremSpheres® are fatigue, nausea and vomiting, abdominal pain, fever, mild to moderate abnormality of liver function tests.

Potential Serious Adverse Events Due to High Radiation

1. **Acute pancreatitis** ---- causes immediate severe abdominal pain. Verify by SPECT imaging of the abdomen and test for serum amylase.
2. **Acute Gastritis** ---- causes abdominal pain. Verify by standard methods to diagnose gastric ulceration.
3. **Acute cholecystitis** ---- causes significant upper abdominal pain and may require cholecystectomy for resolution. Verify by appropriate imaging studies.
4. **Radiation Hepatitis** ---- causes unexplained progressive deterioration of liver function. Verify by transcutaneous core biopsy of the liver.

5. **Radiation Pneumonitis** ---- causes excessive non-productive cough. Verify by x-ray or HR pulmonary CT evidence of pneumonitis.

10. RADIATION SAFETY

The preparation and administration procedure must be regarded as being a potentially serious radiation hazard to the clinical staff. Regulatory and local radiation handling guidelines should be followed concerning set-up, microsphere administration, waste disposal and post-implantation care.

Radiation hygiene principles (ALARA) should be taken into account at all time. In short, this means that dose exposure to clinical staff, nursing staff and unintended dose exposure to the patient should be 'as low as reasonably achievable' by considering the following aspects:

- **TIME** – Minimize the time of exposure
- **DISTANCE** – Increase the distance between the radiation source and body/ body extremities as much as possible
- **SHIELDING** – Take appropriate shielding measures

APPENDIX I – MEASUREMENT OF POST-IMPLANT EXPOSURE

Dose rates were measured at 1.0 m distance from the abdomen of patients (n=3) implanted with a Holmium-166 microspheres activity, associated with a recommended average whole-liver absorbed dose of 60 Gy.

The results of the measurements are presented in Table 2.

Table 2 – Dose rates at 1.0m distance from treated patients at t = 0, 6, 24 and 48h after treatment.

	t = 0h	t = 6h	t = 24h	t = 48h
Dose rate \dot{D} ($\mu\text{Sv/h}$)	31 (18-53)	27 (16-46)	17 (10-29)	9 (5-15)

APPENDIX II – CORRECTION FOR DECAY

The physical half-life of Holmium-166 is 26.8 hours. To calculate the activity of QuiremSpheres® at any moment after the reference time multiply the reference activity with the appropriate decay factor shown in Table 3.

Caution: The reference time must be in the user's local time before correcting for decay.


Table 3 – Decay factors of QuiremSpheres® microspheres

Hours	Decay Factor	Hours	Decay Factor
0.5	0.987	9	0.792
1	0.974	10	0.772
2	0.950	11	0.752
3	0.925	12	0.733
4	0.902	24	0.538
5	0.879	36	0.394
6	0.856	48 (day 2)	0.289
7	0.834	72 (day 3)	0.155
8	0.813	96 (day 4)	0.083

APPENDIX III – EXPLANATION OF SYMBOLS ON PACKAGING AND LABELING

Manufacturer

Temperature limitation

Date of manufacture

Caution, contains radioactive material

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