QuiremScout®

QS-S001

EN	- Instructions for Use	1
DE	- Gebrauchsanweisung	4
ES	- Instrucciones de uso	7
FR	- Mode d'emploi	10
HU	- Használati útmutató	13
IT	- Istruzioni per l'uso	16
LT	- Naudojimo instrukcijos	19
NL	- Gebruiksaanwijzing	22
PL	- Instrukcja użycia	25
PT	- Instruções de utilização	28
SK	- Inštrukcie na používanie	31
SL	- Navodila za uporabo	34
SV	- Bruksanvisning	37
TR	- Kullanım Talimatları	40

QuiremScout®

Quirem Medical B.V.
Zutphenseweg 55
7418 AH Deventer
The Netherlands
www.quirem.com

1. DEVICE DESCRIPTION

QuiremScout® consists of small dose of biocompatible poly(L-lactic acid) (PLLA) microspheres that contain holmium-166. The microspheres have a mean diameter of 25-35 micrometers. Holmium-166 is a beta emitting isotope for therapeutic use. The maximum energy of the beta particles is 1.85 MeV (50.0%) and 1.77MeV (48.7%). The maximum range of the emitted beta particles in tissue is 8.7 mm with a mean of 2.5 mm. In addition, holmium-166 emits primary gamma photons (81 KeV). The half-life is 26.8 hours, which means that more than 90% of the radiation is delivered within the first 4 days following the administration procedure. QuiremScout® is a permanent implant.

QuiremScout® is for single patient use only. Do not re-sterilize the product and do not use the product after the expiry date / time indicated on the product label.

2. PRESENTATION AND STORAGE

QuiremScout® contains a standard dose of approximately 3 million microspheres with an approximate total activity of max. 300 MBq at planned moment of administration. The standard dose is supplied in a single V-vial or divided over multiple (up to three) V-vials. Details of the calibration date, activity at calibration and expiry information are indicated on the package labels. The V-vial(s) comes in a lead container, packed in a type A transport box. The vial and its content should be stored inside its transportation container at room temperature (15-25°C, 59-77°F).

Caution: Upon receiving QuiremScout®, check for package damage. Do not use the product in case the V-vial or the packaging is damaged.

3. INTENDED USE

QuiremScout® is intended for evaluation of lung-shunt, extrahepatic deposition and intrahepatic distribution of intra-arterially injected microspheres.

4. INDICATIONS FOR USE

QuiremScout® is indicated for patients that are eligible for SIRT treatment of liver tumors.

5. PATIENT PREPARATION

The following tests are recommended to be performed prior to administration of QuiremScout®:

- A hepatic angiogram to establish arterial anatomy of the liver.
- Serologic tests of liver function to determine the extent of liver function damage.
- Evaluate the vascular network of the liver using high resolution imaging.

QuiremScout® can be ordered divided over multiple vials to simulate segmental SIRT from multiple injection positions.

Caution: The total administered amount of net activity per patient must not exceed 300 MBq.

6. ACCESORIES

It is highly recommended to use the following accessories for the QuiremScout® administration procedure:

- QuiremSpheres® Delivery Set (QS-D001)
- QuiremSpheres® Customer Kit (QS-C001)

For an extended and illustrated description of the preparation and execution of the QuiremScout® administration procedure, the reader is referred to documentation of the QuiremSpheres® Customer Kit. These instructions should be read and understood in their entirety prior to use.

7. DOSE CALCULATION

The total dose of QuiremScout® (see Section 2) can be ordered in multiple (up to three) vials to be divided amongst multiple injection positions per targeted liver volume.

To calculate the fraction per vial (Fr_i), divide the targeted liver volume for each vial (V_i) by the total targeted liver volume (V_{total}), i.e.:

$$Fr_i = V_i / V_{total}$$

Consequently, all fractions summed should equal 1, i.e.:

$$\sum Fr_i = 1.$$

8. ADMINISTRATION OF QUIREMSCOUT

QuiremScout is administered intra-arterially. Using standard techniques, position a microcatheter into the hepatic artery (common, right, left or one of the other branches to the liver). The microcatheter used for injection must have a minimal internal diameter of at least 0.65 mm.

Slowly inject the microspheres with saline (0.9% NaCl). Do not use an injection rate of more than 5mL saline per minute.

NOTE: When administering QuiremScout®, the catheter tip should be in the same position as planned for the therapeutic SIRT procedure.

9. EVALUATION OF ADMINISTRATION

Assessing the in-vivo QuiremScout® distribution to evaluate lung shunt, extra hepatic deposition and intra-hepatic distribution can be done by gamma scintigraphy or Single Photon Emission Computed Tomography (SPECT). QuiremScout® emits photons with an energy of 81 KeV. To ensure sufficient photon detection, it is recommended to perform imaging on the same day of injection.

10. CONTRAINDICATIONS

QuiremScout® is contraindicated for patients that

- are not eligible for SIRT;
- have uncorrectable abnormal vascular anatomy that would result in significant reflux of hepatic arterial blood to the lungs, stomach, pancreas or bowel;

11. PRECAUTIONS, ADVERSE EVENTS

11.1. General Precautions

- Only clinical staff that are trained under the QuiremScout® training program may order, handle or implant QuiremScout®.
- Safety and effectiveness of this device in pregnant women, nursing mothers or children have not been established.
- This product is radioactive. Local regulations must be followed when handling this device
- This product should be used only by physicians who are qualified by training and experience in the safe use and handling of radionuclides and whose experience and training have been approved by the appropriate government agency authorized to license the use of radionuclides.
- Do not use if the vial or packaging appear damaged.
- Do not use after the expiry date/time indicated on the product label.
- Sterile and single use product. Do not reuse.
- The total administered amount of net activity per patient must not exceed 300 MBq.

11.2. Adverse Events

No adverse events have been reported for QuiremScout® to date.

12. MANAGEMENT OF FORESEEABLE UNWANTED SITUATIONS

Radioactive contamination of staff, equipment or treatment room

In case of (suspected) radioactive contamination of staff, equipment or treatment room:

- Determine the extent of the radioactive contamination by measuring the gamma photon emission of holmium-166 with an appropriate hand-held radiation detector.
- Follow applicable local guidelines on cleaning or quarantining the contaminated surface(s).

Product has expired before time of procedure

Discard the product according to local regulations. Schedule a new procedure and order a new product.

No trained qualified staff is available at time of the procedure

Cancel the procedure and discard the product according to local regulations.

Catheter tip cannot be placed in the desired position

Abort the procedure and discard the product according to local regulations.

13. RADIATION SAFETY

The preparation and administration procedure must be regarded as being a potentially serious radiation hazard to the clinical staff. Regulatory and local radiation hygiene guidelines should be followed concerning set-up, microsphere administration, waste disposal and post-implantation care.

Radiation hygiene principles (ALARA) should be taken into account at all time. In short, this means that dose exposure to clinical staff, nursing staff and unintended dose exposure to the patient should be 'as low as reasonably achievable' by considering the following aspects:

- **TIME** – Minimize the time of exposure.
- **DISTANCE** – Increase the distance between the radiation source and body/ body extremities as much as possible.
- **SHIELDING** – Take appropriate shielding measures.

APPENDIX I – CORRECTION FOR DECAY

The physical half-life of holmium-166 is 26.8 hours. To correct for physical decay of this radionuclide, the fractions that remain at selected intervals after the time of calibration are shown in **Table 1**.

Caution: The time of the calibration indicated on the product label must be converted to local time prior to correcting for decay.

Table 1 – Physical decay table of ¹⁶⁶Ho. Half-life: 26.8 hours

Hours	Fraction remaining	Hours	Fraction remaining
0*	1	9	0.792
0.5	0.987	10	0.772
1	0.974	11	0.752
2	0.950	12	0.733
3	0.925	24	0.538
4	0.902	36	0.394
5	0.879	48 (day 2)	0.289
6	0.856	72 (day 3)	0.155
7	0.834	96 (day 4)	0.083
8	0.813		

*Calibration Time.

APPENDIX II – EXPLANATION OF SYMBOLS

Manufacturer

Date of manufacture

Batch code

Serial number


Use by date

Do not reuse

Consult instructions for use

Temperature limitation

Caution, contains radioactive material

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