TestSpheres

Holmium-166 Microspheres

REF QS-T001

EN - Instructions for Use TestSpheres (QS-T001)	3	LT - Naudojimo instrukcija „TestSpheres“ (QS-T001)	23
CS - Návod k použití TestSpheres (QS-T001)	5	NL - Gebruiksaanwijzing TestSpheres (QS-T001)	25
DA - Brugervejledning TestSpheres (QS-T001)	7	NO - Bruksanvisning TestSpheres (QS-T001)	27
DE - Gebrauchsanweisung TestSpheres (QS-T001)	9	PL - Instrukcja użycia TestSpheres (QS-T001)	29
ES - Instrucciones de uso TestSpheres (QS-T001)	11	PT- Instruções de utilização TestSpheres (QS-T001)	31
FI - Käyttöohjeet TestSpheres (QS-T001)	13	RU – Инструкция по применению TestSpheres (QS-T001)	33
FR - Mode d'emploi TestSpheres (QS-T001)	15	SK - Návod na použitie TestSpheres (QS-T001)	35
HR - Upute za uporabu TestSpheres (QS-T001)	17	SL - Navodila za uporabo TestSpheres (QS-T001)	37
HU - Használati utasítás TestSpheres (QS-T001)	19	SV - Bruksanvisning TestSpheres (QS-T001)	39
IT - Istruzioni per l'uso TestSpheres (QS-T001)	21	TR - Kullanım Talimatları TestSpheres (QS-T001)	41

Manufacturer

Quirem Medical B.V.
Zutphenseweg 55
7418 AH Deventer
The Netherlands

1. INTRODUCTION

This instructions for use is a manual for using TestSpheres in a safe and responsible manner. TestSpheres can be used to verify the response of equipment that is utilized for performing measurements on radioactive Holmium-166 products. This includes performance control of the dose calibrators and measuring detection limits of Single Photon Emission Computed Tomography (SPECT) systems.

2. PRODUCT DESCRIPTION

TestSpheres consists of biocompatible poly-L-lactic acid (PLLA) microspheres containing Holmium-166 suspended in 2 mL fluid. The microspheres have a mean diameter of 30 micrometer. Holmium-166 is a high-energy beta emitting isotope. The maximum energy of the beta particles is 1.85 MeV (50.0%) and 1.77MeV (48.7%). In addition, Holmium-166 emits primary gamma photons (81 KeV, 6.7%). The half-life is 26.8 hours.

TestSpheres is supplied as a standard activity of 2, 4, 6 or 8 GBq at reference time in a single capped V-vial. Details of the calibration date and activity at calibration are indicated on the package labels. The V-vial comes in a lead container, packed in a type A transport box.

3. INTENDED USE

TestSpheres can be used for the following purposes:

- Preparing dose calibrators for measuring Holmium-166 activity
- Performing Holmium-166 saturation measurements of SPECT(-CT) systems

WARNING: TestSpheres is ***not*** for patient use!

4. SAFETY

TestSpheres is radioactive. The preparation, handling and disposal must be regarded as potentially hazardous. Local regulations must be followed and only qualified staff is allowed to handle this product.

Radiation hygiene principles (ALARA) should be considered at all time. In short, this means that dose exposure should be 'as low as reasonably achievable' by considering the following aspects:

- **TIME** – Minimize the time of exposure
- **DISTANCE** – Increase the distance between the radiation source and body/ body extremities as much as possible
- **SHIELDING** – Take appropriate shielding measures

5. INSTRUCTIONS FOR USE

CAUTION: Upon receiving TestSpheres, check for package damage. Do not use the product in case the V-vial or package is damaged.

5.1 Unpacking

- Open the type A transport box
- Take out the lead container
- Remove the sealing tape and take of the lid of the lead container

The V-Vial containing the microspheres can now be taken out for further handling.

WARNING: Never open the V-vial.

CAUTION: If the V-vial has visible damage, be aware of radiation leakage and treat it as such.

CAUTION: When handling the V-vial, always use forceps to maintain distance.

5.2 Handling during measurements

- Take the V-vial out of the lead container using forceps and place it into the device that will be used for measurements.
- After finishing the measurements, put the V-vial back in the lead container using forceps

APPENDIX I – CORRECTION FOR DECAY

The physical half-life of Holmium-166 is 26.8 hours. To calculate the actual activity of TestSpheres, multiply the calibrated activity with the appropriate decay factor.

CAUTION: The time of the calibration mentioned on the product label must be converted to the user's local time.

Table 1 – Decay factors for Holmium-166

Hours	Decay Factor	Hours	Decay Factor
0.5	0.987	9	0.792
1	0.974	10	0.772
2	0.950	11	0.752
3	0.925	12	0.733
4	0.902	24	0.538
5	0.879	36	0.394
6	0.856	48 (day 2)	0.289
7	0.834	72 (day 3)	0.155
8	0.813	96 (day 4)	0.083

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