HolmiumSolution

Holmium-166 Chloride Solution

QS-T002

EN - Instructions for Use HolmiumSolution (QS-T002)	3	LT - Naudojimo instrukcija „HolmiumSolution“ (QS-T002)	23
CS - Návod k použití HolmiumSolution (QS-T002)	5	NL - Gebruiksaanwijzing HolmiumSolution (QS-T002)	25
DA - Brugervejledning HolmiumSolution (QS-T002)	7	NO - Bruksanvisning HolmiumSolution (QS-T002)	27
DE - Gebrauchsanweisung HolmiumSolution (QS-T002)	9	PL - Instrukcja użycia HolmiumSolution (QS-T002)	29
ES - Instrucciones de uso HolmiumSolution (QS-T002)	11	PT- Instruções de utilização HolmiumSolution (QS-T002)	31
FI - Käyttöohjeet HolmiumSolution (QS-T002)	13	RU – Инструкция по применению HolmiumSolution (QS-T002)	33
FR - Mode d'emploi HolmiumSolution (QS-T002)	15	SK - Návod na použitie HolmiumSolution (QS-T002)	35
HR - Upute za uporabu HolmiumSolution (QS-T002)	17	SL - Navodila za uporabo HolmiumSolution (QS-T002)	37
HU - Használati utasítás HolmiumSolution (QS-T002)	19	SV - Bruksanvisning HolmiumSolution (QS-T002)	39
IT - Istruzioni per l'uso HolmiumSolution (QS-T002)	21	TR - Kullanım Talimatları HolmiumSolution (QS-T002)	41

Quirem Medical B.V.
Zutphenseweg 55
7418 AH Deventer
The Netherlands

1. INTRODUCTION

This instructions for use is a manual for using HolmiumSolution in a safe and responsible manner. HolmiumSolution can be used for preparation of phantoms for Single Photon Emission Computed Tomography (SPECT) measurements.

2. PRODUCT DESCRIPTION

HolmiumSolution consists of Holmium-166 chloride dissolved in 2 ml sterile water.. Holmium-166 is a high-energy beta emitting isotope. The maximum energy of the beta particles is 1.85 MeV (50.0%) and 1.77 MeV (48.7%). In addition, Holmium-166 emits primary gamma photons (81 KeV, 6.7%). The half-life is 26.8 hours.

HolmiumSolution is supplied as a standard activity of +/- 600 MBq in a single capped V-vial. Details of the calibration date and activity at calibration are indicated on the package labels. The V-vial comes in a lead container, packed in a type A transport box.

3. INTENDED USE

HolmiumSolution can be used for the following purpose:

- Preparation of phantoms for measuring SPECT(-CT) system related parameters

WARNING: HolmiumSolution is ***not*** for patient use!

4. SAFETY

HolmiumSolution is radioactive. The preparation, handling and disposal must be regarded as potentially hazardous. Local regulations must be followed and only qualified staff is allowed to handle this product.

Radiation hygiene principles (ALARA) should be considered at all time. In short, this means that dose exposure should be 'as low as reasonably achievable' by considering the following aspects:

- **TIME** – Minimize the time of exposure
- **DISTANCE** – Increase the distance between the radiation source and body/ body extremities as much as possible
- **SHIELDING** – Take appropriate shielding measures

5. INSTRUCTIONS FOR USE

CAUTION: Upon receiving HolmiumSolution, check for package damage. Do not use the product in case the V-vial or package is damaged.

5.1 Unpacking

- Open the type A transport box
- Take out the lead container
- Remove the sealing tape and take of the lid of the lead container

The V-vial containing the solution can now be taken out for further handling.

WARNING: Never open the V-vial.

CAUTION: If the V-vial has visible damage, be aware of radiation leakage and treat it as such.

CAUTION: When handling the V-vial, always use forceps to maintain distance.

5.2 Phantom preparation

- 1) Secure the V-vial in an appropriate holder or clip to ensure safe processing.
- 2) Remove the aluminum sealing to open the way to the septum.
- 3) Withdraw the desired amount [ml] of HolmiumSolution from the vial using a syringe or similar device and put it in the solution that will be used to fill the phantom.
- 4) Ensure that the HolmiumSolution is distributed homogenously in the fluid and fill the phantom volume(s) with the fluid.

APPENDIX I – CORRECTION FOR DECAY

The physical half-life of Holmium-166 is 26.8 hours. To calculate the actual activity of HolmiumSolution, multiply the calibrated activity with the appropriate decay factor.

CAUTION: The time of the calibration mentioned on the product label must be converted to the user's local time.

Table 1 – Decay factors for Holmium-166

Hours	Decay Factor	Hours	Decay Factor
0.5	0.987	9	0.792
1	0.974	10	0.772
2	0.950	11	0.752
3	0.925	12	0.733
4	0.902	24	0.538
5	0.879	36	0.394
6	0.856	48 (day 2)	0.289
7	0.834	72 (day 3)	0.155
8	0.813	96 (day 4)	0.083

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